

## PROTOCOL FOR Jnk1 KO MOUSE GENOTYPING

### Procedure

Genotyping of offspring from Jnk1 KO breeding colony is based on PCR.

### PCR primers

5' forward primer (J1-forw) 5' **tcg acc cat gct aag cgc gcc** 3'  
 3' reverse primer (J1-rev) 5' **cta ctt aat aac ggg ggt gga gga tca c**  
 3'  
 3' reverse primer (J1-lacZ) 5' **cgg tgc ggg cct ctt cgc** 3'

### PCR profile – JNKETT

95 °C, 10 min  
 95 °C, 45 s   35 cycles  
 60 °C, 45 s  
 72 °C, 1 min  
 72 °C, 10 min  
 4 °C, ∞

*PCR mix* - the PCR's for KO versus wt are run separate.

10 x PCR Gold buffer (Perkin Elmer)	3.0 μl
MgCl <sub>2</sub> (25 mM)	2.0 μl
dNTPs (10 mM)	0.5 μl
J2-forw (20 μM)	0.5 μl
J2-rev (20 μM)	0.5 μl
J2-lacZ	0.5 μl
AmpliTaq Gold (5 U/μl)	0.2 μl
DNA template (~ 0.5 μg tail DNA)	2.0 μl
ddH <sub>2</sub> O	21.3 μl
	30 μl

### Post-PCR analysis

Load 10 μl of the PCR reaction on a 1,5 % agarose gel.  
 Expected results; two bands – 150 bp (KO; J1-forw + J1-lacZ) and 200 bp  
 (wt; J1-forw + J1-rev)