

D50R genotyping protocol

DNA extraction from tails

- Add 0.5ml tail buffer containing proteinase K to tube with tail.
Tail buffer: 100mM tris pH8
5mM EDTA
0.2% SDS
200mM NaCl
100µg/ml Proteinase K
- Incubate o/n at 50C
- Flick to mix. Centrifuge at 13K (bench top centrifuge) for 20mins
- Remove supernatant to fresh tube
- Add 0.5ml of phenol:chlorophorm:isoamylalcohol
- Flick to mix, centrifuge at 13K for 10min
- Transfer aqueous top layer to new tube
- Add 0.5ml of isopropanol, flick to mix and precipitate DNA
- Centrifuge for 30mins at 13K
- Remove wash (check to see pellet in tube and try not to suck up pellet)
- Add 0.5ml 70% ethanol, flick to wash
- Centrifuge at 13K for 30mins
- Remove supernatant and air dry for 10mins
- Add 40µl of TE
- Keep in fridge use 1 or 2µl in PCR

PCR

The standard PCR reaction mix was set up as shown in Table 2.1. DNA amplification was carried out in a PTC-100 thermal cycler (MJ Research, Inc) using the standard conditions shown in Table 2.2.

Reagent	Volume
10x buffer (Bioline)	5 μ l
Primer 1 (100ng/ μ l)	2.5 μ l
Primer 2 (100ng/ μ l)	2.5 μ l
dNTPs (10mM, Amersham)	1 μ l
taq (Bioline,5U/ μ l)	0.5 μ l
DNA (1 μ g genomic DNA)	1 μ l
H ₂ O	37.5 μ l

Table 2.1. Standard PCR reaction mix.

PCR primers

Primer 1	intron 4aR1	ACTATGAAGTCCAGAGTGCC
Primer 2	intron 2F7	GCATTGAATCCCCAGCATCA

Step	Condition
1	95°C for 5 minutes
2	95°C for 1 minute
3	61°C for 1 minute
4	72°C for 3 minutes
5	Go to step 2 x30
6	72°C for 10 minutes

Table 2.2. Standard PCR conditions

Following PCR, 25 μ l PCR reaction was digested by StuI restriction enzyme digest using standard conditions at 37C for 2 hrs before analysis by agarose gel electrophoresis on a 1.2% gel.

The expected wild type (WT) allele PCR product size was 1470 bp whilst the mutant PCR product sizes were 1133 bp and 334 bp after *StuI* digestion




Figure 5.3. Confirmation of the *Nbr1*^{D50R} knock-in mouse genotype by genomic PCR and sequencing.

A. Schematic diagram showing the location of the PCR primers and *StuI* restriction enzyme sites used for genotyping.

B. Agarose gel electrophoresis separation of PCR product. The PCR product spans from intron 2 to intron 4a and the expected WT PCR product size is 1470bp. The mutant allele contains two introduced *StuI* restriction sites either side of the second *LoxP* site, therefore, on restriction enzyme digest, the mutant PCR product is cut into two sized bands of 1133bp and 334bp. The *LoxP* site is 40bp so cannot be visualised on this gel.

C. Electropherogram showing the correct mutation is present in the D50R knock-in mouse compared with WT (red box).

WT: Wild type; Mut: *Nbr1*^{D50R} mutant.