

Dr Vassilis Aidinis laboratory  
Immunology Research Department  
Alexander Fleming Biomedical Sciences Research Center  
34 Fleming Street,  
Vari, 16672  
Greece  
Tel. +30 210 9656310 (ext. 159)  
Fax +30 210 9653934  
Email: aidinis@fleming.gr

## **Genotyping protocol for lysophosphatidic acid receptor 1 (Lpar1) (MGI ID: 108429) tm1d mice.** (Protocol edited by Evangelos Dioletis on 6-11-2015)

The Lpar1 tm1d mouse was derived by crossing Lpar1 cKO mice [B6-Lpar1<sup>tm1c(EUCOMM)Wtsi</sup>/Flmg] with a Cre driver strain.

The offspring were further mated with C57BL6 mice in order to remove Cre recombinase from the genetic background.


### Genotyping primers

Fwd primer (789): 5'- GGA TGC TAT TCT GGG GAT GA -3' (20 bp, T<sub>m</sub>=58.4°C)  
Rev primer\_1 (790): 5'- ATA CCC AAT GCA GCC AAA AA -3' (20 bp, T<sub>m</sub>=54.3°C)  
Rev primer\_2 (791): 5'- TCA TGG ACA CTT GGA CTA ATG AA -3' (23 bp, T<sub>m</sub>=59.3°C)

The current PCR protocol can discriminate between homo- and heterozygous knockout mice.

A schematic representation of the genotyping strategy along with the primer binding sites and the expected PCR products is shown below:

# Genotyping of Lpar1 tm1d mice



Please note that the expected PCR products for the WT and defloxed (knockout) alleles are 325 bp and 564 bp long respectively.

Typical pipetting scheme for a single reaction (in  $\mu\text{l}$ ):


ddH <sub>2</sub> O	10.4
10x buffer*	2
MgCl <sub>2</sub> (25mM)	1.2
dNTPs (2.5 mM each)	2
Fwd primer (5 $\mu\text{M}$ )	1 (789)
Rev primer 1 (5 $\mu\text{M}$ )	1 (790)
Rev primer 2 (5 $\mu\text{M}$ )	1 (791)
DNA template (50-100 ng/ $\mu\text{l}$ )	1
Taq DNA polymerase	0.4
<b>Total volume (<math>\mu\text{l}</math>)</b>	<b>20</b>

\*10x PCR Buffer: 500 mM KCl, 100mM Tris-HCl (pH 9.0 at 25°C), 1% Triton X-100.

Thermocycler conditions:

95°C for 5 min  
95°C for 30 s  
60.7°C for 30 s } x35  
72°C for 45 s }  
72°C for 5 min

Example of genotyping LPAR1 tm1d mice.



2% agarose gel electrophoresis in 1x TBE  
35 min at 120V (6.7V/cm).  
M: Marker (*Pst*I-cut lambda DNA)

## ***Pst*I-cut lambda DNA**




Figure taken from: <http://www.taq-dna.com/phage-lambda-dna-psti-digest-ready-to-use-144.html>