

**Name of Mouse model or mutation:****Prps1-R214W-EM1-B6****Prps1-R214W-EM2-B6****Description:**

Point mutation generated using CRISPR/Cas9 reagents.

**Type of mutation:**

Point mutation: R214W

**Sequence details****Prps1 WT:**

```
GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTCTAGACTCACAAGTC
TTCCATGATGGAAATTGGGATAATCTGTACTCCCTGCACTAGCTATAAGATCTCTCTATCTTCTG
CTTTTAGTTAGCAAGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTAC
CCCATCTATACTATGGACTGTGCAAGTTTGAATATTATATTCTATACTGCTGTTGGTATTACAGAA
GCATCAATCCACAAAATCCAGTTTGTATTGTTTGTGGGATGTTAATGATAGTTGCTTTCTTCTA
GCTCAGGGATCTGTCTTACTTTCAAGGAGTATTCTCCTACAGTAAAAACATTCGGTTATTTTGTGTA
GTGGGCACAAGCGATGGAACCTGGGGCCAGCTGTCTACCACTAAATGACACTCCCAGCTTTGGAAG
CCGCTTGCCTTGGAATTTTGTCTTTCCTAACTTTGAACATCTATTTTACAGAGTGACCTCCATTGCAG
ACCGGCTGAATGTGGATTTTGTCTTGTATTACAAAGGAACGGAAGAAGGCCAATGAAGTGGACCGCA
TGGTACTTGTGGGAGATGTGAAGGACCGTGTGGCTATACTTGTGGATGACATGGCTGACACTTGCG
GTACAATCTGTCACGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAATAAAAAGCTAACTACTCTT
TACTGATTCTTCTGAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTCTGATCATG
TTAACTTTTAATAACATATAGATTAGTAAGTATTCTCATTTCAGATTTTTTTAAAGATTTATTTTATA
TCTAAATACACTGTAGCTGTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTTG
TGAGCCACCATGTGGTTGCTGGGATTTGTACTCAGGACCTTTGGAAGAGCAGCCAGTGCTCTTAACT
GCTGAGCCATCTCTCCAGCCCTCATTTCAGATTTTTAACTTTTATAAACAGACTGCTTTTACATCCA
GGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTT
AGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTACATTTCTGCCTGGAAGCTGTGTC
AAGACAGTATATTGATTGCATATTGCTAGACCAGGGAAAGATGAAA
```


**Prps1-R214W-EM1-B6 and Prps1-R214W-EM2-B6:**

```
GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTCTAGACTCACAAGTC
TTCCATGATGGAAATTGGGATAATCTGTACTCCCTGCACTAGCTATAAGATCTCTCTATCTTCTG
CTTTTAGTTAGCAAGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTAC
CCCATCTATACTATGGACTGTGCAAGTTTGAATATTATATTCTATACTGCTGTTGGTATTACAGAA
GCATCAATCCACAAAATCCAGTTTGTATTGTTTGTGGGATGTTAATGATAGTTGCTTTCTTCTA
```

GCTCAGGGATCTGTCTTACTTTCAAGGAGTATTCTCCTACAGTAAAAACATTCCGTTATTTTGTGTAA  
 GTGGGCACAAGCGATGGAACCTGGGGCCAGCTGTCTACCACTAAATGACACTCCCAGCTTTGGAAG  
 CCGCTTGCCTTGGATTTTGTCTTTCCTAACTTTGAACATCTATTTTACAGAGTGACCTCCATTGCAG  
 ACCGGCTGAATGTGGATTTTGCTTTGATTCACAAGGAACGGAAGAAGGCCAATGAAGTGGACCGCA  
 TGGTACTTGTGGGAGATGTGAAGGACTTGGGTGGCTATACTTGT**CG**ATGACATGGCTGACACTTGCG  
 GTACAATCTGTCACGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAATAAAAAGCTAAACTACTCTT  
 TACTGATTCTTTCTGAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTCTGATCATG  
 TTAACCTTTAATAACATATAGATTAGTAAGTATTCTCATTTCAGATTTTTTTAAAGATTTATTTTTATA  
 TCTAAATACACTGTAGCTGTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTTG  
 TGAGCCACCATGTGGTTGCTGGGATTTGTAICTCAGGACCTTTGGAAGAGCAGCCAGTGTCTTAACT  
 GCTGAGCCATCTCTCCAGCCCTCATTTCAGATTTTTAACTTTTATAAACAGACACTGCTTTTACATCCA  
 GGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTT  
 AGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTACATTTCTGCCTGGAAGCTGTGTC  
 AAGACAGTATATTGATTGCATATTGCTAGACCAGGGAAAGATGAAA

Nucleotide changes highlighted in **red and underlined = nominated change**, silent changes highlighted in red only.

**Prps1-R214W-EM1-B6 and Prps1-R214W-EM2-B6 Heterozygous F1 animal sequence trace:**



### Prps1-R214W-EM1-B6 and Prps1-R214W-EM2-B6 Nucleotide Alignment:

```

*      20      *      40      *      60      *      80      *      100     *      120     *      140     *
Prps1_WT      : GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTTAGACTCACAAGCTTCCATGATGGAATTTGGGATAATCTGTACACTCCCTGCCTAGCTATAAGATCTCTCTATCTTCTGCTTTTAGTTAGCA
Prps1_R214W_EM1 : GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTTAGACTCACAAGCTTCCATGATGGAATTTGGGATAATCTGTACACTCCCTGCCTAGCTATAAGATCTCTCTATCTTCTGCTTTTAGTTAGCA
Prps1_R214W_EM2 : GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTTAGACTCACAAGCTTCCATGATGGAATTTGGGATAATCTGTACACTCCCTGCCTAGCTATAAGATCTCTCTATCTTCTGCTTTTAGTTAGCA
GAGGGTAGGCTTTACTCCACTCCTGCCACCTTCAGGTTCACTTGTCACTAATTTAGACTCACAAGCTTCCATGATGGAATTTGGGATAATCTGTACACTCCCTGCCTAGCTATAAGATCTCTCTATCTTCTGCTTTTAGTTAGCA

      160      *      180      *      200      *      220      *      240      *      260      *      280      *      300
Prps1_WT      : AGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTACCCCATCTATATACTATGGACTGTGCAAGTTTGAATATTATATTCCTATACTGCTGTTGGTATTACAGAAGCATCAATCCACAAAATCCAGTTTGT
Prps1_R214W_EM1 : AGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTACCCCATCTATATACTATGGACTGTGCAAGTTTGAATATTATATTCCTATACTGCTGTTGGTATTACAGAAGCATCAATCCACAAAATCCAGTTTGT
Prps1_R214W_EM2 : AGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTACCCCATCTATATACTATGGACTGTGCAAGTTTGAATATTATATTCCTATACTGCTGTTGGTATTACAGAAGCATCAATCCACAAAATCCAGTTTGT
AGCTAATATACCTTAGTATATGCTGCTCCAGGACTGGTGAAGGGTATGCCTACCCCATCTATATACTATGGACTGTGCAAGTTTGAATATTATATTCCTATACTGCTGTTGGTATTACAGAAGCATCAATCCACAAAATCCAGTTTGT

*      320      *      340      *      360      *      380      *      400      *      420      *      440      *
Prps1_WT      : ATTGTTTGTTTTGTGGGATGTTAATGATAGTTGCTTCTCTAGCTCAGGGATCTGCTTACTTTCAAGGAGTATTCCTCCACAGTAAAAACATTCCTGTTATTTTGTGTAAGTGGGCACAAGCGATGGAACCTGGGGCCAGCTGCTTACC
Prps1_R214W_EM1 : ATTGTTTGTTTTGTGGGATGTTAATGATAGTTGCTTCTCTAGCTCAGGGATCTGCTTACTTTCAAGGAGTATTCCTCCACAGTAAAAACATTCCTGTTATTTTGTGTAAGTGGGCACAAGCGATGGAACCTGGGGCCAGCTGCTTACC
Prps1_R214W_EM2 : ATTGTTTGTTTTGTGGGATGTTAATGATAGTTGCTTCTCTAGCTCAGGGATCTGCTTACTTTCAAGGAGTATTCCTCCACAGTAAAAACATTCCTGTTATTTTGTGTAAGTGGGCACAAGCGATGGAACCTGGGGCCAGCTGCTTACC
ATTGTTTGTTTTGTGGGATGTTAATGATAGTTGCTTCTCTAGCTCAGGGATCTGCTTACTTTCAAGGAGTATTCCTCCACAGTAAAAACATTCCTGTTATTTTGTGTAAGTGGGCACAAGCGATGGAACCTGGGGCCAGCTGCTTACC

      460      *      480      *      500      *      520      *      540      *      560      *      580      *      600
Prps1_WT      : ACTAAATGACACTCCCAGCTTTGGAAGCCGCTTTGCCTTGGATTTTGTCTTTCCCTAATTTGAAACATCTATTTTACAGAGTGACCTCCATTGCAGACCCGGCTGAATGTGGATTTTGCCTTTGATTACAAAGGAACGGAAGGCCAATG
Prps1_R214W_EM1 : ACTAAATGACACTCCCAGCTTTGGAAGCCGCTTTGCCTTGGATTTTGTCTTTCCCTAATTTGAAACATCTATTTTACAGAGTGACCTCCATTGCAGACCCGGCTGAATGTGGATTTTGCCTTTGATTACAAAGGAACGGAAGGCCAATG
Prps1_R214W_EM2 : ACTAAATGACACTCCCAGCTTTGGAAGCCGCTTTGCCTTGGATTTTGTCTTTCCCTAATTTGAAACATCTATTTTACAGAGTGACCTCCATTGCAGACCCGGCTGAATGTGGATTTTGCCTTTGATTACAAAGGAACGGAAGGCCAATG
ACTAAATGACACTCCCAGCTTTGGAAGCCGCTTTGCCTTGGATTTTGTCTTTCCCTAATTTGAAACATCTATTTTACAGAGTGACCTCCATTGCAGACCCGGCTGAATGTGGATTTTGCCTTTGATTACAAAGGAACGGAAGGCCAATG

*      620      *      640      *      660      *      680      *      700      *      720      *      740      *
Prps1_WT      : AAGTGGACCCGATGGTACTTGTGGGAGATGTGAAGGACGGTGGGCTATACTTGTGATGACATGGCTGACACTTGCCTGACAAATCTGTACAGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAAATAAAAGCTAAACTACTCTTTACTG
Prps1_R214W_EM1 : AAGTGGACCCGATGGTACTTGTGGGAGATGTGAAGGACGGTGGGCTATACTTGTGATGACATGGCTGACACTTGCCTGACAAATCTGTACAGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAAATAAAAGCTAAACTACTCTTTACTG
Prps1_R214W_EM2 : AAGTGGACCCGATGGTACTTGTGGGAGATGTGAAGGACGGTGGGCTATACTTGTGATGACATGGCTGACACTTGCCTGACAAATCTGTACAGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAAATAAAAGCTAAACTACTCTTTACTG
AAGTGGACCCGATGGTACTTGTGGGAGATGTGAAGGACGGTGGGCTATACTTGTGATGACATGGCTGACACTTGCCTGACAAATCTGTACAGCCGCTGACAAGTGAGTACAGGCTAGTTAGAAAATAAAAGCTAAACTACTCTTTACTG

      760      *      780      *      800      *      820      *      840      *      860      *      880      *      900
Prps1_WT      : ATTCTTTCTGAAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTTCTGATCATGTTAACTTTTAATAACATATAGATTAGTAAGTATTCCTCATTTCAGATTTTTTTAAAGATTTATTTTTATATCTAAATACACTGTAGCT
Prps1_R214W_EM1 : ATTCTTTCTGAAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTTCTGATCATGTTAACTTTTAATAACATATAGATTAGTAAGTATTCCTCATTTCAGATTTTTTTAAAGATTTATTTTTATATCTAAATACACTGTAGCT
Prps1_R214W_EM2 : ATTCTTTCTGAAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTTCTGATCATGTTAACTTTTAATAACATATAGATTAGTAAGTATTCCTCATTTCAGATTTTTTTAAAGATTTATTTTTATATCTAAATACACTGTAGCT
ATTCTTTCTGAAAACACAACCTGAAGGGTAAGTAGAGTAGTAGCTTTGGGATATTTCTGATCATGTTAACTTTTAATAACATATAGATTAGTAAGTATTCCTCATTTCAGATTTTTTTAAAGATTTATTTTTATATCTAAATACACTGTAGCT

*      920      *      940      *      960      *      980      *      1000     *      1020     *      1040     *
Prps1_WT      : GTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTGTGAGCCACCATGTGGTGTCTGGGATTTGTAAGTACAGGACCTTTGGAAGAGCAGCCAGTGTCTTAACTGCTGAGCCATCTCTCCAGCCCTCATTTCAGATTT
Prps1_R214W_EM1 : GTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTGTGAGCCACCATGTGGTGTCTGGGATTTGTAAGTACAGGACCTTTGGAAGAGCAGCCAGTGTCTTAACTGCTGAGCCATCTCTCCAGCCCTCATTTCAGATTT
Prps1_R214W_EM2 : GTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTGTGAGCCACCATGTGGTGTCTGGGATTTGTAAGTACAGGACCTTTGGAAGAGCAGCCAGTGTCTTAACTGCTGAGCCATCTCTCCAGCCCTCATTTCAGATTT
GTCTTCAGATGCACTAGAAGAGGGCGTCAGATCTCATTATAGCTGGTGTGAGCCACCATGTGGTGTCTGGGATTTGTAAGTACAGGACCTTTGGAAGAGCAGCCAGTGTCTTAACTGCTGAGCCATCTCTCCAGCCCTCATTTCAGATTT

      1060     *      1080     *      1100     *      1120     *      1140     *      1160     *      1180     *      1200
Prps1_WT      : TTAACCTTTTATAAACAAGACACTGCTTTTACATCCAGGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTTAGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTAC
Prps1_R214W_EM1 : TTAACCTTTTATAAACAAGACACTGCTTTTACATCCAGGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTTAGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTAC
Prps1_R214W_EM2 : TTAACCTTTTATAAACAAGACACTGCTTTTACATCCAGGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTTAGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTAC
TTAACCTTTTATAAACAAGACACTGCTTTTACATCCAGGATATCCTTAATGAGTTATCAAAAGTTAAATCCATCACATACACTTTACCTTTTAAAAGCCCTTACTTAGTAACACAATCTACCATAGACTGGGTTGTCTCCCTTCATGATTAC

*      1220     *      1240     *      1260
Prps1_WT      : ATTTCTGCCTGGAAGCTGTGTCAAGACAGTATATTTGATTGCATATTTGCTAGACCAGGGAAGATGAAA
Prps1_R214W_EM1 : ATTTCTGCCTGGAAGCTGTGTCAAGACAGTATATTTGATTGCATATTTGCTAGACCAGGGAAGATGAAA
Prps1_R214W_EM2 : ATTTCTGCCTGGAAGCTGTGTCAAGACAGTATATTTGATTGCATATTTGCTAGACCAGGGAAGATGAAA
ATTTCTGCCTGGAAGCTGTGTCAAGACAGTATATTTGATTGCATATTTGCTAGACCAGGGAAGATGAAA
```

**Prps1-R214W-EM1-B6 and Prps1-R214W-EM2-B6 Predicted Protein Alignment:**

```
                *      20      *      40      *
Prps1_WT       : VTSIADRLNVDFALIHKERKKANEVDRMVLVGDVKDRAVAILVDDMADTCGTICHAADK*
Prps1_R214W_EM1 : VTSIADRLNVDFALIHKERKKANEVDRMVLVGDVKDRAVAILVDDMADTCGTICHAADK*
Prps1_R214W_EM2 : VTSIADRLNVDFALIHKERKKANEVDRMVLVGDVKDRAVAILVDDMADTCGTICHAADK*
                VTSIADRLNVDFALIHKERKKANEVDRMVLVGDVKDRAVAILVDDMADTCGTICHAADK
```

### QC strategy employed at Harwell to check the edited allele:

Genomic DNA was extracted from ear clip biopsies and amplified in a PCR reaction using the following conditions/primer sequences:

Geno_PRPS1_R214W_F4	GAGGGTAGGCTTTACTCCAC
Geno_PRPS1_R214W_R4	TTTCATCTTCCCTGGTCTAGCAA
Taq Polymerase used	Roche Expand Long Range DNTPack
Annealing Temperature (°C)	64
Elongation time (min)	1
WT product size (bp)	1268
Mutant product size (bp)	1268
Notes	Sequencing primers:  Geno_Prps1_R214W_F3 AGGACTGGTGAAGGGTATGC  Geno_Prps1_R214W_R3 GGGAGACAACCCAGTCTATGGTA

All amplicons were sent for Sanger sequencing to check for integration of the donor oligo sequence at the target site. F1 sequences should be heterozygous unless on sex chromosome.

As the sgRNA used had a potential off-target site with two or fewer mismatches, then PCR amplification and Sanger sequencing of this site was performed to confirm its integrity using the following conditions/primer sequences:

Geno_Prps1_OT3_F1	CTGTAGCAGGTGCGGATCAT
Geno_Prps1_OT3_R1	GCAATGGAGCATGAGCTTTGT
Taq Polymerase used	Platinum SuperFi PCR Master Mix
Annealing Temperature (°C)	58
Elongation time (min)	1
WT product size (bp)	927
Mutant product size (bp)	n/a

Notes	
-------	--

Copy counting of the donor sequence was carried out by ddPCR at the F1 stage to confirm donor oligos were inserted once on target into the genome. The following Taqman assay was used to copy count the donor sequence compared against a VIC-labelled reference assay for Dot1l:

Assay name	Prps1-R214W-donor-UNIV1
Forward Primer (5'-3')	ACTTGTGGGAGATGTGAAGGAC
Reverse Primer (5'-3')	GCTTTTATTTCTAACTAGCCTGTA CTAC
Probe (5'-3')	ATGACATGGCTGACACTTGCGGTA
Label	FAM-BHQ1

The ddPCR assay recognises sequence common to both the WT Prps1 and the R214W mutant. Therefore, WT controls and correctly targeted F1 R214W heterozygote animals will call at 2 copies.

Assay name	Prps1-R214W-donor-MUT1
Forward Primer (5'-3')	GTGGGAGATGTGAAGGACTGG
Reverse Primer (5'-3')	GCTTTTATTTCTAACTAGCCTGTA CTCA
Probe (5'-3')	TGGCTATACTTGTCGATGACATGGC
Label	FAM-BHQ1

The ddPCR assay is specific to the Prps1 R214W mutation. Therefore, WT controls will call at 0 copies and correctly targeted F1 R214W heterozygote animals will call at 1 copy.

Reference Assay Name	Dot1l
Forward primer (5'-3')	GCCCCAGCACGACCATT
Reverse primer (5'-3')	TAGTTGGCATCCTTATGCTTCATC
Probe (5'-3')	CCCAACAGGCCTGGATTCTCAATGC
Label	VIC

VIC-labelled reference assay for Dot1l gene.



# PRPS1-R214W-EM1-B6

# PRPS1-R214W-EM2-B6

## Allele Description

This is a CRISPR/Cas9 induced mutation creating point mutation R214W in exon ENSMUSE00000384936 of *Prps1*. The stock was generated at MRC Harwell via microinjection of CRISPR/Cas9 reagents into 1-cell stage embryos.

## qPCR Copy Counting Genotyping Strategy

The genotyping strategy presented here has been optimized for reagents and conditions used by the Genotyping Core at MRC Harwell. To genotype animals, we recommend researchers validate the assay independently. PCR cycling temperature and times may require additional optimization based on the specific genotyping reagents used.

Samples are genotyped using qPCR copy counting with both a wild type and a mutant assay against a known reference assay (*Dot1l* on chromosome 10; 2 copies present). Samples for this line are genotyped using the following primers and probe:

- Universal probe and Universal primer designed 5' of the point mutation.
- Wild type (WT) specific primer binding to the WT base mutated in the mutant allele.
- Mutant specific primer binding to the R214W point mutation.

For X-linked genes that have been targeted, the following results would be expected:

Genotype of the Modified allele	WT Assay	Mutant Assay
Wildtype male	1	0
Hemizygous mutant male	0	1
Wildtype female	2	0
Heterozygous female	1	1
Homozygous mutant female	0	2



# PRPS1-R214W-EM1-B6

# PRPS1-R214W-EM2-B6

## Prps1-R214W-WT1 assay (FAM labelled)

ACATCTATTTTACAGAGTGACCTCCATTG**CAGACCGGCTGAATGTGGATT**TTTGCTTTGATTCACAAGGAACGG  
**AAGA**AGGCCAATGAAGTGGACCGCATGGTACTTGTGGGAGATGTGAAGGAC**Cg**t**GTGGCTATACTTGTg**GAT  
GACATGGCTGACACTTGC~~GG~~TACAATCTGT**CAC**CGCGCTGACAAGT**GAGTACAGG**CTAGTTAGAAATAAA

Lower case letters denote bases changed in the mutant allele.  
Probe sequence is in bold and shaded grey.  
Primer sequences are in bold and underlined.

Oligo Name	5' label	Sequence 5' → 3'	3' label	Oligo Type
Prps1-R214W-UNI_F	n/a	<u>CAGACCGGCTGAATGTGGATT</u>	n/a	Universal Forward
Prps1-R214W-UNI_PROBE	FAM	<b>TGCTTTGATT</b> CACAAGGAACGG <b>AAGA</b>	BHQ	Universal Probe
Prps1-R214W-WT_R	n/a	<u>CATCCACAAGTATAGCCACACG</u>	n/a	Wild type Reverse

## Prps1-R214W-MUT1 assay (FAM labelled)

ACATCTATTTTACAGAGTGACCTCCATTG**CAGACCGGCTGAATGTGGATT**TTTGCTTTGATTCACAAGGAACGG  
**AAGA**AGGCCAATGAAGTGGACCGCATGGTACTTGTGGGAGATGTGAAGGAC**tGg****GTGGCTATACTTGTc**GAT  
GACATGGCTGACACTTGC~~GG~~TACAATCTGT**CAC**CGCGCTGACAAGT**GAGTACAGG**CTAGTTAGAAATAAA

Lower case letters denote bases changed in the mutant allele.  
Probe sequence is in bold and shaded grey.  
Primer sequences are in bold and underlined.

Oligo Name	5' label	Sequence 5' → 3'	3' label	Oligo Type
Prps1-R214W-UNI_F	n/a	<u>CAGACCGGCTGAATGTGGATT</u>	n/a	Universal Forward
Prps1-R214W-UNI_PROBE	FAM	<b>TGCTTTGATT</b> CACAAGGAACGG <b>AAGA</b>	BHQ	Universal Probe
Prps1-R214W-MUT_R	n/a	<u>CATCGACAAGTATAGCCACCCA</u>	n/a	Mutant Reverse





# PRPS1-R214W-EM1-B6

# PRPS1-R214W-EM2-B6

## Dot1l internal control (VIC labelled)

CTGATGGGTGTGGGCAGATCCTACAGAGTCCCATTGGCCACCATGTGTGCTACGCCTGAAATAAAGCCTT**GCC**  
**CCAGCACGACCATT**CAGGG**CCAGCTCTCAAGTCG**ACTGTAA**GATGAAGCATAAGGATGCCAACTA**CTAACA  
 GAAAACGACTAGAGGGGAAAAGAACAAGGAAACAGAAGACGCAGCACTCCGGCTTCCCTGGGTTGGCCAGT  
 CACCCATGA

Oligo Name	5' label	Sequence 5' → 3'	3' label	Oligo Type
Dot1l_Forward	n/a	<b><u>GCCCCAGCACGACCATT</u></b>	n/a	WT Forward
Dot1l_Probe	VIC	<b>CCAGCTCTCAAGTCG</b>	BHQ	WT Probe
Dot1l_Reverse	n/a	<b><u>TAGTTGGCATCCTTATGCTTCATC</u></b>	n/a	WT Reverse

Probe sequence is in bold and shaded grey

Primer sequences are in bold and underlined

## DNA extraction method

DNA is extracted from ear clips using Applied Biosystems Taqman Sample-to-SNP Kit and qPCR run using 1:10 dilution from the crude preparation.

## qPCR master mix 1X

Applied Biosystems GTX Taqman master mix	5 µl
Dot1l_Forward (20 µM)	0.225 µl
Dot1l_Reverse (20 µM)	0.225 µl
Dot1l_Probe (5 µM)	0.2 µl
FAM Assay (probe 5 µM & primers 15 µM each)	0.3 µl
ddH2O	1.55 µl
DNA (1:10 dilution of ABI Sample-to-SNP prep)	2.5 µl

Each sample is ran in technical duplicate. Seven WT and/or mutant controls are also included in duplicate along with non-template controls.

## qPCR cycling conditions

qPCR instrument: Applied Biosystems 7500/7900 or ThermoFisher QuantStudio 7

95°C for 20 sec  
 Then 40 cycles of;  
 95°C for 3 sec  
 60°C for 30 sec





# PRPS1-R214W-EM1-B6 PRPS1-R214W-EM2-B6

Version No. 1  
Date: 18.02.2020  
Created/Updated by: AC