

## PCR FOR CRE (common)

### PCR REACTION:

10X PCR buffer	<b>2ul</b>
MgCl <sub>2</sub> 25 mM	Depends on Taq polymerase ( <b>usually 2ul</b> )
dNTP'S 2,5 mM	1ul
Primer A	5 pmol per reaction ( <b>1ul from 5pmole/ul stock</b> )
Primer B	5 pmol per reaction ( <b>1ul from 5pmole/ul stock</b> )
DNA	1ul
Taq	Variable ( <b>0.4ul</b> )
dd.H <sub>2</sub> O	Up to volume (20λ) <b>11.6ul</b>

**Primer A (sense) → 30:** 5'-att-acc-ggt-cga-tgc-aac-gag-t-3'

**Position in relating gene:** (bacteriophage P1 cre recombinase )68-80\*

**Primer B (antisense) → 31:** 5'-cag-gta-tct-ctg-acc-aga-gtc-a-3'

**Position in relating gene:**(bacteriophage P1 cre recombinase )852-873\*

(\*starting point: atg of cre recombinase)

### PCR CONDITIONS:

1. 94°C for 5 min.
2. 94°C for 1 min.
3. 57°C for 1 min.
4. 72°C for 1 min.
5. Steps 2-4 for 5 more cycles.
6. 94°C for 10 sec.
7. 55°C for 40 sec.
8. 72°C for 1:30 min.
9. Steps 6-8 for 25 more cycles.
10. 72°C for 10 min.
11. 16°C for 1 min.
12. End

**PCR PRODUCT:** Transgene band: ~800bp.

**Samples:**

### PHOTO EXAMPLE:

