



Genotyping protocol

Acap3

IR00002865 / E136

(ICS internal reference)

This report has been prepared by: **Valérie Rousseau**
33 (0)3 88 65 56 55
genotyping@igbmc.fr

This report has been validated by: **Sylvie Jacquot, PhD, Head of Genotyping Service**
33 (0)3 88 65 57 44
genotyping@igbmc.fr

The first version of this report was finalized the: 05 Dec 2012

For any question, please contact:

Institut Clinique de la Souris - ICS - Mouse Clinical Institute
1 rue Laurent Fries, BP 10142
67404 Illkirch Cedex, France
Email: genotyping@igbmc.fr
Web site: <http://www-mci.u-strasbg.fr/>

TABLE OF CONTENTS

Table of contents2

1. Genotyping protocol and data2

 1.1. Genotyping strategy2

 1.2. PCR protocol4

 1.3. Picture of genotyping with various alleles4

2. Cre and Flp genotyping method5

 2.1. Cre and Flp genotyping5


 2.2. PCR Protocol6

1. Genotyping protocol and data

This section describes the condition used at the Mouse Clinical Institute (ICS) to genotype your **Acap3** Constitutive Knockout / Conditional Knockout (KO-cKO) project.

1.1. Genotyping strategy

The map below describes the position of the primers used for genotyping for each possible allele.



Sequence of primers used for genotyping:

| Position | Primers | Sequence |
|----------|---------|-----------------------------|
| Ef | 5221 | ATGAGCTCCAAGGAGACCCGATG |
| Ef2 | 5219 | TGCCTTTGGCTTCCGGGATT |
| Er | 5223 | TGTTGGTGCCCATGATAGTCTTCTCTG |
| Kr | 3210 | CCTGTCCCTCTCACCTTCTACC |
| L3f | 5222 | TTAAGGAATGTCTGCAGAGGTTGGAG |
| L3r | 5220 | CGCTGACGTCCAGAAATAAGTTGTCCT |
| Lxr | 5086 | GAAGTTATCATTAATTGCGTTGCGCC |

PCR fragments expected size (bp):

| Region analyzed | Primers used | Position on the primer (see the map above) | Targeted allele (KO allele) (L3) | cKO allele (L2) | KO allele (L-) | WildType allele (WT) |
|--|--------------|---|----------------------------------|-----------------|----------------|----------------------|
| 5' part of the selection marker | 5221-3210 | Ef / Kr | 461 | --- | --- | --- |
| Presence of the distal loxP | 5222-5220 | L3f / L3r | 360 | 360 | --- | 413 |
| Distal loxP specific PCR | 5222-5086 | L3f / Lxr | 245 | 245 | --- | --- |
| Excision of the selection marker | 5219-5223 | Ef2 / Er | 7452* | 548 | --- | 368 |
| Excision of the floxed exon(s), i.e. knock out | 5219-5220 | Ef2 / L3r | 8899* | 1995* | 491** | 1868** |

*: this PCR product will not be observed using our PCR genotyping conditions (see description below)

** : this PCR is only verified if mice are generated

---: no Amplicon should be obtained

1.2. PCR protocol

This section describes the composition of the mix and cycling conditions used for genotyping.

| Reagents: | Volume: |
|--------------------------------|-------------|
| - FastStart PCR Master (Roche) | 7.5µl |
| - DNA (50ng/µl) | 1.5µl |
| - 5' primer (100 µM) | 0.06µl |
| - 3' primer (100 µM) | 0.06µl |
| - Sterile H ₂ O | up to 15 µl |

Cycling conditions:

| Temp | Time | #Cycles |
|------|------|---------|
| 95°C | 4min | 1 |
| 94°C | 30s | 34 |
| 62°C | 30s | |
| 72°C | 1min | |
| 72°C | 7min | 1 |
| 20°C | 5min | 1 |

NB: These PCR conditions have been optimized for high-throughput genotyping. Adaptation to small-scale may be required.


2. Cre and Flp genotyping method

The protocol used to segregate the cre and/or flp transgene is indicated below.

Detection of cre transgene and flp transgene is done using a multiplex assay: primer pairs were designed for each gene and for a positive control (Cpxm1 gene).

2.1. Cre and Flp genotyping

Schematic representation of the genotyping strategy



Sequence of primers used for genotyping:

| Primers | Sequence |
|---------|---------------------------|
| 4029 | ACTGGGATCTTCGAACTCTTTGGAC |
| 4030 | GATGTTGGGGCACTGCTCATTACC |
| 4045 | CCATCTGCCACCAGCCAG |
| 4046 | TCGCCATCTTCCAGCAGG |
| 4195 | TCTTTAGCGCAAGGGGTAGGATCG |
| 4196 | GTCCTGGCCACGGCAGAAGC |

PCR fragments expected size (bp):

| Primer pair | 4045-4046 | 4195-4196 | 4029-4030 |
|-----------------|------------------------------|------------------------------|--------------------|
| Region analyzed | Middle part of Cre transgene | Middle part of Flp transgene | Cpxm1 control gene |
| Control gene | / | / | 397 |
| Tg allele | 281 | 328 | / |

2.2. PCR Protocol

This section describes the composition of the mix and cycling conditions used for genotyping.

| Reagents | Volume |
|------------------------------|-------------|
| FastStart PCR Master (Roche) | 7.5µl |
| DNA (50ng/µl) | 1.5µl |
| 5' primer (100 µM) | 0.05µl |
| 3' primer (100 µM) | 0.05µl |
| Sterile H ₂ O | up to 15 µl |

Cycling conditions are identical to those described in chapter 1.2