



Genotyping protocol

Dcaf7

/ P5838

(ICS internal reference)

This report has been prepared by: **David Moulaert**
genotyping@igbmc.fr

This report has been validated by: **Sylvie Jacquot, PhD, Head of Genotyping Service**
33 (0)3 88 65 57 44
genotyping@igbmc.fr

The first version of this report was generated the: 28 Apr 2016

For any question, please contact:

Institut Clinique de la Souris - ICS - Mouse Clinical Institute
1 rue Laurent Fries, BP 10142
67404 Illkirch Cedex, France
Email: genotyping@igbmc.fr
Web site: <http://www-mci.u-strasbg.fr/>

TABLE OF CONTENTS

Table of contents2

1. Genotyping protocol and data2

 1.1. Genotyping strategy2

 1.2. PCR protocol4

2. qPCR Genotyping protocol and data5

 2.1. Genotyping strategy5

 2.2. qPCR protocol6


3. Cre and Flp genotyping method7

1. Genotyping protocol and data

This section describes the condition used at the Mouse Clinical Institute (ICS) to genotype your **Dcaf7** Constitutive Knockout / Conditional Knockout (KO-cKO x Cre) project.

1.1. Genotyping strategy

The map below describes the position of the primers used for genotyping for each possible allele.



Sequence of primers used for genotyping:

Position	Primers	Sequence
Ef	8900	CCTGTGTGCACAGTGAGCTTGAGG
Kr	3209	CCAACAGCTTCCCCACAACGG
L3r	8899	CTCCCCTCGCTTTCCTGATAAACAG
Lxr	3255	ACTGATGGCGAGCTCAGACCATAAC
Rif	5966	GCACATGGCTGAATATCGACGGT

²: for a selected position, a second primer was designed

PCR fragments expected size (bp):

Region analyzed	Primers used	Position on the primer (see the map above)	Targeted allele (HR)	KI allele	WildType allele
5' part of the selection marker	8900-3209	Ef / Kr	416	---	---
Cre total excision	5966-3255	Rif / Lxr	5040*	471	---
Cre total excision	5966-8899	Rif / L3r	5104*	535	---

*: this PCR product will not be observed using our PCR genotyping conditions (see description below)

** : this PCR is only verified if mice are generated

---: no Amplicon should be obtained

1.2. PCR protocol

This section describes the composition of the mix and cycling conditions used for genotyping.

Reagents:	Volume:
- FastStart PCR Master (Roche)	7.5µl
- DNA (50ng/µl)	1.5µl
- 5' primer (100 µM)	0.06µl
- 3' primer (100 µM)	0.06µl
- Sterile H ₂ O	up to 15 µl

Cycling conditions:

Temp	Time	#Cycles
95°C	4min	1
94°C	30s	34
62°C	30s	
72°C	1min	
72°C	7min	1
20°C	5min	1


NB: These PCR conditions have been optimized for high-throughput genotyping. Adaptation to small-scale may be required.

2. qPCR Genotyping protocol and data

This section describes the condition used at the Mouse Clinical Institute (ICS) to genotype your **Bscl2** Constitutive Knockout / Conditional Knockout (KO-cKO x Cre) project.

2.1. Genotyping strategy

The map below describes the position of the primers used for genotyping for each possible allele.



Sequence of primers used for genotyping:

Position	Primers	Sequence
Rf	7443	CTCGCCACTTCAACATCAAC
Rr	7445	TTATCAGCCGGAAAACCTACC
Mf	Neof1	TGAATGAACTGCAGGACGAG
Mr	Neor1	TCCCGCTTCAGTGACAAC

2.2. qPCR protocol

Reagents:	Volume:
- EvaGreen (biorad)	3,5µl
- DNA (10ng/µl)	3µl
- Forward primer (100µM)	0,06µl
- Reverse primer (100µM)	0,06µl
- Sterile H2O	up to 7µl

Cycling conditions:

Temp	Time	#Cycles
95°C	10min	1
95°C	5s	
62°C	10s	34
95°C	15min	
Melting curve analysis		
65°C -> 95°C		

Follow manufacturer's protocol for programming the data acquisition of dsDNA product.

NB: These PCR conditions have been optimized for high-throughput genotyping. Adaptation to small-scale may be required.

3. Cre and Flp genotyping method

You will find the genotyping protocol in the publication:

[Highly-efficient, fluorescent, locus directed cre and FlpO deleter mice on a pure C57BL/6N genetic background.](#)

Birling MC, Dierich A, Jacquot S, Hérault Y, Pavlovic G.

Genesis. 2012 Jun;50(6):482-9. doi: 10.1002/dvg.20826. Epub 2012 Mar 20.