

STRAIN NAME : **ThBcl2.1**

REAGENTS AND SOLUTIONS:

- Kapa Mouse Genotyping kit (Sigma)
- Primers 30µM (Eurofins)
- Mineral Oil (Fisher Scientific)
- DNA (extracted with Kapa Mouse Genotyping kit)

REACTION MIX :

MIX	Quantity / well
Water	7,40 µl
2X Kapa Buffer	10,00 µl
ThBcl2F	0,30 µl
ThBcl2R	0,30 µl
DNA	2,00 µl
Final Quantity	20,00 µl

PCR PROGRAM:

3 min	95°C	
15 sec	94°C	
15 sec	60°C	35 cycles
30 sec	72°C	
5 min	72°C	

MIGRATION :

Migration on 1.5% agarose gel.

PRIMERS SEQUENCES

- ThBcl2F 5' CCTTCTTTGAGTTCGGTGGGG 3'
- ThBcl2R 5' GCCAAACTGAGCAGAGTCTTCAGAG 3'

EXPECTED AMPLIFICATIONS

Tg allele : 210 bp

Perfection date : 21/11/2018

Update :