

PCR AMH-cre

1- PCR mix

2,50 µl 10X buffer(Q BIOgene ; 1,5 mM MgCl₂)
 0,50 µl dNTP (20 mM total)
 0,15 µl primer AMH1 (200 ng/µl)
 0,15 µl primer AMH2 (200 ng/µl)
 0,20 µl Taq polymerase (Q BIOgene ; 5U/µl)
 20,50 µl distilled water

2- PCR reaction :

24,0 µl mix
 1,0 µl DNA (200ng/µl)

3- PCR program :

5 min	94°C	
1 min	94°C	
1 min	55°C	35 cycles
2 min	72°C	
5 min	72°C	

4- Primers:

- Primer AMH1 5' TCCT GGA AAA TGC TTC TGT CCG 3'
 - Primer AMH2 5' CAG GGT GTT ATA AGC AAT CCC 3'

5- Amplifications:

transgene + : 400 bp