

Wellcome Trust Genome Campus Hinxton Cambridge CB10 1SA T 01223 834244 F 01223 494919 www.sanger.ac.uk

Gene: Sdc3

Colony prefix: DANG

Allele: Sdc3 em1(IMPC)Wtsi

## Allele type: Crispr/Cas9 mediated Point mutation

## Allele information: T329

Further information about the allele can be found on the 'International Mouse Phenotyping Consortium' (IMPC) web site at <u>http://www.mousephenotype.org/</u>

| 33460<br>Sd | 33500 | 33520 | 33540 | 33560 | 33580 | 33600 | 33620 | 33640 | 33660 | 33680 | 33700 | 33720 | 33740 | 33760 | 33780 | 33800<br>SdR |
|-------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------------|
|             |       |       |       |       |       | C     | G>ta  |       |       |       |       |       |       |       |       |              |

# Mouse QC information

| Loss of WT Allele<br>(LOA) qPCR | na | Mutation Sequence<br>confirmed  | pass |
|---------------------------------|----|---------------------------------|------|
| Mutant Specific SR-<br>PCR      | na | Off-target analysis<br>complete | na   |

## Guide RNAs and mutant oligos used in initial experiment

| Sequence   | Chr | Chr Start | Chr End   |
|--|-----|-----------|-----------|
| CCGGACACGGCCAATGAGGTGGT  | 4   | 130821473 | 130821495 |
| CCAGTAAGTGGGGGGGCCCAGCGGGGACTTTGAGCTTC<br>AAGAAGAGACCACGCAGCCGGACATAGCCAATGAGGT<br>GGTGGCTGTGGAAGGAGCCGCGGCCAAGCCGTCACCT<br>CCACTGGGGACA | 4   | 130821419 | 130821541 |

This technical data sheet and information ("Datasheet") is supplied by Genome Research Limited ("GRL").



## Mutant allele sequence:

# Genotyping by end-point PCR

## PCRs primer pairs and expected size bands

| Assay Type   | Assay      | Forward Primer | Reverse Primer | Expected Size Band (bp) |
|--------------|------------|----------------|----------------|-------------------------|
| Standard PCR | Screening* | Sdc3_PM_WT_F   | Sdc3_PM_WT_R   | 341                     |

\*The screening PCR flanks the SNP region and can be used for sequence verification of the allele. The PCR will not distinguish wild type from mutant mice, however, as a product will be amplified in all cases.

# We recommend that mice are sequence-verified with the screening primers to confirm the genotyping qPCR results when establishing the colony, in case of any cross-talk between the assays.

#### **Primer sequences**

| Primer Name  | Primer Sequence (5' > 3') |
|--------------|---------------------------|
| Sdc3_PM_WT_F | GCGTACCTAACCCCAAGTCA      |
| Sdc3_PM_WT_R | GGAGCCTGCCCTCACCTAC       |

#### **Reaction setup**

| Reagent                   | μl   |
|---------------------------|------|
| DNA (~50-100 ng)          | 1    |
| 10x Buffer                | 2    |
| MgCl2 (50 mM)             | 0.6  |
| Platinum Taq (Invitrogen) | 0.2  |
| dNTPs (100 mM)            | 0.2  |
| Primer 1 (10 μM)          | 0.4  |
| Primer 2 (10 μM)          | 0.4  |
| ddH20                     | 15.2 |
| Total                     | 20   |

This technical data sheet and information ("Datasheet") is supplied by Genome Research Limited ("GRL").



Wellcome Trust Genome Campus Hinxton Cambridge CB10 1SA T 01223 834244 F 01223 494919 www.sanger.ac.uk

### **Amplification conditions**

| Step | Conditions            | Time     |
|------|-----------------------|----------|
| 1    | 94°C                  | 5 min    |
| 2    | 94°C                  | 30 sec   |
| 3    | 58°C                  | 30 sec   |
| 4    | 72°C                  | 1:30 sec |
| 5    | Go to '2' + 34 cycles | -        |
| 6    | 72°C                  | 5 min    |
| 7    | 12°C                  | Forever  |

# Genotyping by SNP qPCR

## Primers for LoA qPCR assay

| Gene | Source       | Forward Primer Seq. | Reverse Primer Seq. | Probe Primer Seq.    |
|------|--------------|---------------------|---------------------|----------------------|
| Sdc3 | Life         | GCCAGAGGTGCCAGTAA   | CCGCGGCTCCTTCCA     | [VIC]CGGACACGGCCAAT  |
|      | Technologies | GTG                 |                     | [FAM]CCGGACATAGCCAAT |

Reactions are performed in a 10µl volume using an Applied Biosystems 7900HT Fast Real-Time PCR System or Applied Biosystems Viia7 with DNA prepared using the Sample-to-SNP<sup>™</sup> kit (Applied Biosystems) from mouse ear biopsies. GTXpress<sup>™</sup> buffer is also used (Applied Biosystems).

| Reagent              | μl   |
|----------------------|------|
| 2x GTXpressTM buffer | 5    |
| 40x target assay     | 0.25 |
| ddH2O                | 3.75 |
| DNA                  | 1    |

#### **Amplification conditions**

| Step      | Conditions     | Time   |
|-----------|----------------|--------|
| Pre-read  | 60°C           | 30 sec |
| 1         | 95°C           | 20 sec |
| 2         | 95°C           | 10 sec |
| 3         | 60°C           | 30 sec |
| 4         | Go to '2' + 34 | -      |
| Post-read | 60°C           | 30 sec |

This technical data sheet and information ("Datasheet") is supplied by Genome Research Limited ("GRL").



Wellcome Trust Genome Campus Hinxton Cambridge CB10 1SA T 01223 834244 F 01223 494919 www.sanger.ac.uk

## Links to information and frequently asked questions

MGP mouse phenotype data: http://www.mousephenotype.org

How the "critical" exon is decided: http://www.i-dcc.org/kb/entry/102/

# **Relevant publications**

White, J.K., Gerdin, A.-K., Karp, N.A., Ryder, E., Buljan, M., Bussell, J.N., Salisbury, J., Clare, S., Ingham, N.J., Podrini, C., et al. (2013). Genome-wide Generation and Systematic Phenotyping of Knockout Mice Reveals New Roles for Many Genes. Cell 154, 452–464.

Mali P, Yang L, Esvelt KM, et al (2013) RNA-guided human genome engineering via Cas9. Science 339:823–6. doi: 10.1126/science.1232033

Jinek M, Chylinski K, Fonfara I, et al (2012) A programmable dual-RNA-guided DNA endonuclease in adaptive bacterial immunity. Science 337:816–21. doi: 10.1126/science.1225829

Cong L, Ran FA, Cox D, et al (2013) Multiplex genome engineering using CRISPR/Cas systems. Science 339:819–23. doi: 10.1126/science.1231143

Singh P, Schimenti JC, Bolcun-Filas E (2014) A Mouse Geneticist's Practical Guide to CRISPR Applications. Genetics genetics.114.169771–. doi: 10.1534/genetics.114.169771

Brandl C, Ortiz O, Röttig B, et al (2015) Creation of targeted genomic deletions using TALEN or CRISPR/Cas nuclease pairs in one-cell mouse embryos. FEBS Open Bio 5:26–35. doi: 10.1016/j.fob.2014.11.009

Zhou J, Wang J, Shen B, et al (2014) Dual sgRNAs facilitate CRISPR/Cas9 mediated mouse genome targeting. FEBS J. doi: 10.1111/febs.12735

Kraft K, Geuer S, Will AJ, et al (2015) Deletions, Inversions, Duplications: Engineering of Structural Variants using CRISPR/Cas in Mice. Cell Rep. doi: 10.1016/j.celrep.2015.01.016

Shen B, Zhang J, Wu H, et al (2013) Generation of gene-modified mice via Cas9/RNA-mediated gene targeting. Cell Res 23:720–3. doi: 10.1038/cr.2013.46

Wang H, Yang H, Shivalila CS, et al (2013) One-step generation of mice carrying mutations in multiple genes by CRISPR/Cas-mediated genome engineering. Cell 153:910–8. doi: 10.1016/j.cell.2013.04.025

Yang H, Wang H, Shivalila CS, et al (2013) One-Step Generation of Mice Carrying Reporter and Conditional Alleles by CRISPR/Cas-Mediated Genome Engineering. Cell 154:1370–1379. doi: 10.1016/j.cell.2013.08.022

This technical data sheet and information ("Datasheet") is supplied by Genome Research Limited ("GRL").